



**SYMMETRIC
OCHRATOXIN WINE**

LATERAL FLOW TEST KIT

for the quantitative determination of Ochratoxin in Wine, must and grape juice

ProGnosis Biotech S.A. is ISO 9001:2015 certified by TÜV Hellas (TÜV NORD).

Use only the current version of Product Data Sheet enclosed with the kit.

Symmetric OCHRATOXIN Wine, S6224/S6248, is a Lateral Flow Test kit for the quantitative determination of Ochratoxin A in wine.

This kit contains all reagents required for 24 or 48 reactions.

Matrices:

Type I: Wine, must, grape juice

- Test time (reaction time after samples and reagents preparation): 5min
- Range: 0 - 5ppb
- Shelf life: 12 months
- Storage: 2-8°C

Specifications

- The LOD of the method is 0.2ppb OTA.
- The LOQ of the method is 0.35ppb OTA.
- Cross-reactivity: The cross-reaction of the anti-Ochratoxin antibody with Ochratoxin A and B is 100 and <0.1% respectively.

1. Description

Symmetric OCHRATOXIN Wine is an innovative Lateral Flow test, utilizing state-of-the-art features for the quantitative detection of Ochratoxin A (OTA) in wine and must. This Lateral Flow test does not require any extraction step. Samples above 2ppb of OTA can also be detected by visual observation, without the use of S-flow reader.

2. General Information

Ochratoxins are a group of mycotoxins produced by some *Aspergillus* species (mainly *A. ochraceus*, but also by 33% of *A. niger* industrial strains) and some *Penicillium* species, especially *P. verrucosum* and *P. carbonarius*. Ochratoxin A (OTA) is the most prevalent and relevant fungal toxin of this group, while ochratoxins B and C are of lesser importance. OTA is a potent nephrotoxin and causes both acute and chronic effects in the kidneys of all mammalian species tested. It is also genotoxic (damages DNA) and teratogenic (damages the fetus) and is considered a probable carcinogen, causing renal carcinoma and other cancers in a number of animal species. Most controlling government agencies worldwide have regulations regarding the amount of aflatoxins allowable in human and animal foodstuffs. Accurate and rapid determination of the presence of OTA in commodities is of paramount importance.

3. Principle of the Method

The quantitative lateral flow test is based on the immunochromatography assay principles. The wells of the microtiter strips contain OTA specific antibodies conjugated to colloidal gold. Diluted sample is added into the well. A dipstick with two capture lines, test and control, is dipped into the well. The suspended mixture starts flowing vertically on the dipstick and passes through the two lines. While running, OTA (if it is present) binds to the antibodies. A valid test should always have the upper control line red. If the sample is free of OTA, a color development occurs at the test line, indicating the absence of OTA in the sample. On the contrary, the presence of OTA in the sample will cause a reduced colored signal at the test line. The test line color intensity is indirectly proportionate to the concentration of OTA present in the samples. By utilizing S-Flow software and the symmetric quantification technology, OTA is accurately quantified.

4. Reagents Provided

Symmetric OCHRATOXIN Wine kit contains sufficient reagents and materials for 24/48 reactions.

Reagents (Store at 2-8°C)	Quantity for 24 reactions	Quantity for 48 reactions
Pots each with 1 strip of 8 reagent microwells and 8 dipsticks	3	6
Sample Diluent Tubes	24	48

5. Materials required but not provided

- 200 or 300µl adjustable micropipettes (single or multi channel) with disposable tips
- S-Flow software along with matching scanner device

6. Storage Instructions

Store kit components between 2 - 8°C. Do not freeze any components provided. Reseal the unused strips in the storing tube together with the desiccant bag provided. The expiry date of the kit and reagents is stated on their labels and no quality guarantee is accepted after the expiration date. The expiry of the kit components can only be guaranteed if the components are stored properly and the reagent is not contaminated due to prior handling. Do not interchange individual components between kits of different lot numbers.

7. Safety and Precautions for use

All reagents should be brought to room temperature (21 - 25°C) before use (at least half an hour) and covered when not in use. Use a clean disposable plastic pipette tip for each reagent, to avoid cross contamination.

8. Sample preparation


Add 200µl of sample into the Sample Diluent Tube provided and mix well.

9. Method Procedure

1. Before opening the reagents, take the kit out of the fridge and wait until the temperature of the reagents reaches the ambient temperature.
2. Open one plastic pot and take out as many test strips and microwells as samples to be tested.
3. The pot with dipsticks should **always be well closed** after reagents have been taken out.
4. Dispense **200µl of diluted sample** into the microwell and pipette **up and down 4 times** to completely mix the lyophilized gold particles in the sample, while avoiding bubbles. The sample should turn into a **uniform pink color**. In case of more than 2 samples, an 8 channel multipipette should be used.
5. Place the appropriate number of sticks into microwells immediately.
6. When the 5 minutes are over, take the dipsticks out of the microwells.
7. Remove the white cotton sample-pad of the stick. Touch the stick with your hand from the colorful pad and remove the white pad with your hands. Do not use a paper, towel or any other material.
8. Place the stick inside the plastic holder in order to be scanned. In case of S-Flow scanner, the sticks must be facing up. In case of EPSON scanner, the **sticks must be facing down (inverted)** and the colored side must be facing the orange sticker.
9. Use S-flow software to quantify results within 10 minutes after the end of analysis. The software will use a Lot specific curve to calculate the results (ppb) according to the matrix sample type.

10. Visual Interpretation

When the test time is completed (5 min), the test strip can also be visually read and interpreted according to the following figure.



Samples containing 2ppb or more of Ochratoxin, can be easily detected via optical observation of the sticks. While test line in negative samples is denser (more reddish) than control line, in samples containing 2ppb of Ochratoxin, test line is more than 50% lighter than control line. In samples with 4ppb or more of Ochratoxin, test line can barely be detected.

11. Performance Evaluation

11.1 Reference Materials


Several reference materials are being used for the evaluation of each product of ProGnosis Biotech S.A. in the context of Quality Control performed by Quality Control Department. Please request a validation report, including the results, at info@prognosis-biotech.com.

11.2 Proficiency Tests

All products participate frequently in Proficiency Tests. For more information, visit the individual product page in our website: www.prognosis-biotech.com

12. Method Summary

Total method time: 5 minutes



VERSION N2

CAT.NUMBER: S6224/S6248

STORAGE: 2-8°C



LATERAL FLOW TEST KIT

for the quantitative determination of Ochratoxin in Wine, must and grape juice

ProGnosis Biotech S.A. is ISO 9001:2015 certified by TÜV Hellas (TÜV NORD).

Use only the current version of Product Data Sheet enclosed with the kit.

Symmetric OCHRATOXIN Wine, S6224/S6248, is a Lateral Flow Test kit for the quantitative determination of Ochratoxin A in wine.

This kit contains all reagents required for 24 or 48 reactions.

Matrices:

Type I: Wine, must, grape juice

- **Sample preparation:** extraction
- **Test time** (reaction time after samples and reagents preparation): 5min
- **Range:** 0 - 5ppb
- **Shelf life:** 12 months
- **Storage:** 2-8°C

All immune assays supplied by ProGnosis Biotech S.A., are warranted to meet or exceed our published specification when used under normal conditions in your laboratory. If the product fails during the stated period, a replacement product will be issued.

ProGnosis Biotech S.A. makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. ProGnosis Biotech S.A. shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product. This method is considered to be a screening method, before a legal action, samples detected as positives must be confirmed with a confirmation method. This product is meant to be used only For Research or Manufacturing use and by qualified technicians.



www.prognosis-biotech.com
e: info@prognosis-biotech.com
t: +30 2410 623922
f: +30 700 700 6262
Farsalon 153 | 41335 Larissa, Greece

